

Nesting Cup

Dose controlled nesting material



- Nesting cups are manufactured using high quality fine kraft paper.
- Packed minimum 1000, preformed cups per box.
- Preformed Nesting cups weigh 8 to 10gm each.
- You know the fixed costs per cage with each use.
- Developed to mimic natural nesting materials, such as thin stemmed grass.
- Ensures that mice can fulfil their natural instincts.
- Helps mice with their thermoregulation and reduces stress.
- Autoclavable.
- Analysis available



For more information:

Tel: +44 (0)1293 827940

Email: sales@lbs-biotech.com

LBS (Serving Biotechnology) Ltd. PO Box 431, RH6 0UW

www.lbs-biotech.com