

A12

Animal marker Pens

The Animal Marker was developed for easy, long lasting identification of laboratory animals.

It is intended for use on the fur, not on the skin.

- Nontoxic and devoid of any hazardous materials.
- A single marker can be used to mark 300-500 animals (depending upon animal size and/or applied area).
- Once applied, colours remain for approximately 12 weeks.
- Broad and Fine tip Markers available.



Available in:

Black
Red
Blue
Green
Orange
Yellow



For more information:

Tel: +44 (0)1293 827940

Email: sales@lbs-biotech.com

LBS (Serving Biotechnology) Ltd. PO Box 431, RH6 0UW

www.lbs-biotech.com