

Pure - Water Gel



Certain animal and protocol conditions create a requirement for “special dietary needs” products to ensure animals’ survival and to encourage reliable research outcomes

Bio-Serv was the first to create a product to meet these special needs in 2003, with the introduction of Nutra-Gel – a highly palatable, nutritionally complete food and water gelled diet

This has now been expanded to include a water gel and an electrolyte gel

These products can be easily fed to animals, by placing the cup or pouch at cage level. This allows the animal to remain healthy and hydrated even through impaired conditions

Pure-Water Gel is a clear gelled water source that can be used as a complete water replacement or supplement for mobility impaired research animals



Gelled, clear water source for lab animals



Provides water without the spills



Convenient cup allows for total replacement or supplemental hydration



Economical and easy to use requiring no refrigeration or special storage



Suitable for a variety of research species



Gamma Irradiated



12 month shelf life

Pyramid Feeder

This well designed stable feeder is ideal for use with critical care diets and specifically with the Pure-WaterGel, ElectroGel and NutraGel products

- ✓ Certified (Contaminant Screened)
- Made of durable pure virgin nylon
- ✗ Use with any powdered, liquid, soft or gelled diet
- ✚ Provides easy access to critical care diets for sick or debilitated animals
- ✓ Autoclavable and cagewasher safe



All Gels can be warmed up and it will turn back to a liquid, you can then add any medication or drugs to it, and it will reset back to a gel!