

Product of the Month...

June 2017



CLIPPER-HOME Exclusive European Distributor



CLIPPER-HOME



- ❑ Made of high temp polycarbonate
- ❑ Offers rodents shelter and security
- ❑ Transparency makes routine cage checks easy and fast
- ❑ Helps reduce boredom and aggressive behaviour
- ❑ Increases physical activity
- ❑ Increases usable floor space, clips to top of cage.
- ❑ Autoclavable
- ❑ Cagewasher safe
- ❑ Certificate of Analysis available

Dimensions

92mm Long x 62mm Width x 57mm depth



For more information:

Tel: +44 (0)1293 827940

Email: sales@lbs-biotech.com

LBS (Serving Biotechnology) Ltd. PO Box 431, RH6 0UW

www.lbs-biotech.com