

Product of the Month...

May 2019



ALPHA-twist™

The Premiere Rodent Enrichment Bedding and Nesting Material

ALPHA-twist™ is comprised of twisted multi-layered paper that is chopped into short lengths. Once in a cage, rodents unfurl and manipulate them into fluffy nests.



Features and Benefits:

- **Maximum Environmental Enrichment:** Species specific behaviors are nurtured and enhanced through the unrolling of the material for use in nest building. One roll of ALPHA-twist™ can be unrolled into a 5 inch strand of paper!
- **Customized Dosage:** Technicians can custom dose each cage with an exact amount of material to satisfy a particular cage or housing arrangement. Other cotton based nesting materials are pre-cut into individual doses which limit flexibility and cost reduction procedures.
- **Economical:** One bag of ALPHA-twist™ can provide enrichment for over 1400 standard cages. The cost per cage for exceptional enrichment is less than traditional cotton based nesting enrichment items.
- **Safer to Use:** ALPHA-twist is safer to use due to its low dust content and the absence of long cotton fiber strands eliminating strangulation issues.

ALPHA-twist™ is also offered pre-blended in our PLUS line of bedding materials:

ALPHA-dri® + PLUS

&

Shepherd's Cob + PLUS



For more information:

Tel: +44 (0)1293 827940

Email: sales@lbs-biotech.com

LBS (Serving Biotechnology) Ltd. PO Box 431, RH6 0UW

www.lbs-biotech.com