

Airocide Air Purification

airocide®

DEVELOPED BY **NASA**

Reduce the risk of infection by destroying airborne bacteria, viruses, mold, and more with Airocide HD.



NASA Developed
Technology



Completely GREEN,
No Emissions

FDA


FDA Class II Listed
Medical Device



Kills 99.99% of
Airborne Pathogens

This incredible and diverse NASA Developed air purification technology can be used to help reduce Respiratory Tract Infections and stop COVID 19 in its tracks!

Utilizing a NASA Developed Photo Catalytic Oxidation Process (PCO), anything organic in the air – from Viruses, Bacteria, Mold/Fungus, Allergens, Ethylene Gas, and other Volatile Organic Compounds (VOC's) get mineralized in AiroCide HD's Bio-Conversion reactor bed with trace elements of CO₂ and water vapor. This technology is a FDA Class II listed Medical Device – and is completely safe and user friendly. The technology does not use nor produce any type of ozone or other chemical/gas – all the kill happens inside the reactor bed with no Emissions.



The Airocide aids in the elimination of animal allergens and smells by destroying:

	VOCs Ethylenes, Formaldehydes, Phenols
	Viruses Influenza, Swine flu, Pneumonia, SARS-CoV-2
	Spores, Mold Cladosporium, Penicillium Citrinum, Aspergillus Niger
	Bacteria Tuberculosis, MRSA, Legionella
	Mycotoxins Aflatoxin, Ochratoxin, Citrinin

Airocides are the ideal choice for holding rooms, surgeries, media rooms, laboratories and cage wash areas.

Please get in touch to see what an Airocide can do for you and your department!



For more information:

Tel: +44 (0)1293 827940

Email: sales@lbs-biotech.com

LBS (Serving Biotechnology) Ltd. PO Box 431, RH6 0UW

www.lbs-biotech.com